



Product Manual

CYTO-ID[®] Autophagy Detection Kit 2.0

Catalog No. ENZ-KIT175

ENZ-KIT175-0200

200 flow cytometry assays

250 fluorescence microscopy assays

3 x 96-well microplate assays

ENZ-KIT175-0050

50 flow cytometry assays

60 fluorescence microscopy assays

1 x 96-well microplate assays



Product Manual

USE FOR RESEARCH PURPOSES ONLY

Unless otherwise specified expressly on the packaging, all products sold hereunder are intended for and may be used for research purposes only and may not be used for food, drug, cosmetic or household use or for the diagnosis or treatment of human beings. Purchase does not include any right or license to use, develop or otherwise exploit these products commercially. Any commercial use, development or exploitation of these products or development using these products without the express written authorization of Enzo Life Sciences, Inc. is strictly prohibited. Buyer assumes all risk and liability for the use and/or results obtained by the use of the products covered by this invoice whether used singularly or in combination with other products.

The CYTO-ID® Autophagy Detection Kit 2.0 is a member of the CELLESTIAL® product line, reagents and assay kits comprising fluorescent molecular probes that have been extensively benchmarked for live cell analysis applications. CELLESTIAL® reagents and kits are optimal for use in demanding cell analysis applications involving confocal microscopy, flow cytometry, microplate readers and HCS/HTS, where consistency and reproducibility are required.

LIMITED WARRANTY; DISCLAIMER OF WARRANTIES

These products are offered under a limited warranty. The products are guaranteed to meet all appropriate specifications described in the package insert at the time of shipment. Enzo Life Sciences' sole obligation is to replace the product to the extent of the purchasing price. All claims must be made to Enzo Life Sciences, Inc., within five (5) days of receipt of order. THIS WARRANTY IS EXPRESSLY IN LIEU OF ANY OTHER WARRANTIES OR LIABILITIES, EXPRESS OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, AND NON- INFRINGEMENT OF THE PATENT OR OTHER INTELLECTUAL PROPERTY RIGHTS OF OTHERS, AND ALL SUCH WARRANTIES (AND ANY OTHER WARRANTIES IMPLIED BY LAW) ARE EXPRESSLY DISCLAIMED.

TRADEMARKS AND PATENTS

Several Enzo Life Sciences products and product applications are covered by US and foreign patents and patents pending. Enzo, CELLESTIAL and CYTO-ID are registered trademarks of Enzo Life Sciences, Inc.

**FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

TABLE OF CONTENTS

INTRODUCTION	2
AUTOPHAGY PATHWAY SCHEMATIC	4
SAFETY WARNINGS AND PRECAUTIONS	5
REAGENTS PROVIDED AND STORAGE	5
ADDITIONAL MATERIALS REQUIRED	6
METHODS AND PROCEDURES	6
A. Reagent Preparations	7
1) Positive Controls.....	7
2) 1X Assay Buffer.....	8
3) CYTO-ID® Green Detection Reagent 2.....	8
B. Live Cell Analysis by Fluorescence/Confocal Microscopy (Adherent Cells).....	9
C. Live Cell Analysis by Fluorescence/Confocal Microscopy (Suspension Cells).....	10
D. Live Cell Analysis by Flow Cytometry.....	11
E. Live Cell Analysis by Fluorescence Microplate Reader ..	12
APPENDICES	15
Fluorescence Channel Selection	15
Typical Results	15
Fluorescence Microscopy	17
Flow Cytometry	18
Fluorescence Microplate	19
Autophagy Treatments	20
REFERENCES	21
TROUBLESHOOTING GUIDE	23
CONTACT INFORMATION	26



Please read entire booklet before proceeding with the assay.



Carefully note the handling and storage conditions of each kit component.



Please contact Enzo Life Sciences Technical Support if necessary.

INTRODUCTION

When subjected to certain hostile conditions that threaten survival, such as when extracellular nutrients are limiting, eukaryotic cells employ a lysosome-mediated intracellular bulk degradation pathway for digesting their own cellular contents by a process referred to as autophagy. Various cytoplasmic constituents, including organelles and long-lived proteins, are sequestered into double-membraned autophagosomes, which subsequently fuse with lysosomes where their contents are degraded. Under physiological conditions, autophagy plays a variety of important roles including maintenance of the amino acid pool during starvation, damaged protein and organelle turnover, prevention of neurodegeneration, tumor suppression, cellular differentiation, clearance of intracellular microbes and regulation of innate and adaptive immunity. Autophagy is considered to be a dynamic, multi-step process which can be regulated at several steps, in both a positive and negative manner. Autophagic activity is typically low under basal conditions, but can be markedly up-regulated, both in cultured cells and intact organisms, by a variety of physiological stimuli such as nutrient starvation, hypoxia, energy depletion, endoplasmic reticulum stress, elevated temperature, high density growth conditions, hormonal stimulation, pharmacological agent treatment, innate immune signaling, and in diseases such as viral, bacterial or parasitic infections as well as various protein aggregopathies (e.g., Alzheimer's, Huntington's and Parkinson's disease), heart disease and acute pancreatitis. Autophagy can be suppressed in certain other diseases, including particular types of cancers, neuro-degenerative disorders, infectious diseases, and inflammatory bowel disorders. A reduction in autophagic function is also considered a characteristic of the aging process.

A conventional fluorescent probe, monodansylcadaverine (MDC), has served as a useful fluorescent marker for lysosomal/ autophagic vacuoles. However, it is known to generate high background and weak fluorescent signal. Enzo Life Sciences' CYTO-ID® Autophagy Detection Kit 2.0 has been optimized for detection of autophagy in live cells by fluorescence microscopy, flow cytometry and fluorescence microplate assay. The assay provides a rapid, specific and quantitative approach for monitoring autophagic activity at the cellular level. The 488 nm-excitable green fluorescent detection reagent supplied in the CYTO-ID® Autophagy Detection Kit 2.0 becomes brightly fluorescent in vesicles produced during autophagy and has been validated under a wide range of conditions known to modulate autophagy pathways.



Product Manual

The conventional way of monitoring autophagic activity is to measure the increased numbers of autophagosomes in cells responding to induction. However, the autophagosome formation is an intermediate stage in the whole dynamic autophagy process. The accumulation of autophagosomes can represent either an increased generation of autophagosomes or the blocked conversion to autolysosomes. To distinguish these two sources of autophagosomes, monitoring autophagic flux provides a meaningful way, in which the balance of the autophagosome generation and clearance will be measured¹⁻⁴.

Rapamycin and starvation are well known inducers of autophagy. Chloroquine is a lysosomal inhibitor. Rapamycin and Chloroquine are included as positive controls in the kit. A nuclear counterstain is provided in the kit as well to highlight cellular nuclei. This live cell analysis kit provides a convenient approach for the analysis of the regulation of autophagic activity at the cellular level.

The CYTO-ID[®] Autophagy Detection Kit 2.0 contains an improved version of the dye found in CYTO-ID[®] Autophagy Detection Kit (ENZ-51031). The mechanism of action remains the same, while the dye has been optimized to increase brightness and photostability.




Figure 1. Schematic depiction of autophagy.⁽⁴⁾ Cytosolic material is sequestered by an expanding membrane sac, the phagophore, resulting in the formation of a double-membrane vesicle, an autophagosome. The outer membrane of the autophagosome subsequently fuses with the lysosome, and the internal material is degraded in the autolysosome. Various regulators of autophagy are also depicted in the diagram.

SAFETY WARNINGS AND PRECAUTIONS

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.



Handle with care

- Reagents should be treated as possible mutagens and should be handled with care and disposed of properly.
- Observe good laboratory practices. Gloves, lab coat, and protective eyewear should always be worn. Never pipet by mouth. Do not eat, drink or smoke in the laboratory areas. All blood components and biological materials should be treated as potentially hazardous and handled as such. They should be disposed of in accordance with established safety procedures.



Protect from light

- To avoid photobleaching, perform all manipulations in low light environments or protected from light by other means.
- The CYTO-ID[®] Green Detection Reagent 2 and the Autophagy Inducer (Rapamycin) contain DMSO which is readily absorbed through the skin. DMSO is harmful if ingested or absorbed through the skin and may cause irritation to the eyes. Observe appropriate precautions when handling these reagents.

REAGENTS PROVIDED AND STORAGE



Avoid freeze / thaw cycles

All reagents are shipped on dry ice. Upon receipt, the kit should be stored at -20°C, protected from light. When stored properly, these reagents are stable for at least twelve months. **Avoid repeated freezing and thawing.**

Reagents provided in the **ENZ-KIT175-0050** kit are sufficient for approximately 50 flow cytometry, 60 fluorescence microscopy or 1 x 96-well microplate assays.

Reagents provided in the **ENZ-KIT175-0200** kit are sufficient for approximately 200 flow cytometry, 250 fluorescence microscopy or 3 x 96-well microplate assays.

Reagent	ENZ-KIT175-0050	ENZ-KIT175-0200
CYTO-ID [®] Green Detection Reagent 2	13 µL	50 µL
Hoechst 33342 Nuclear Stain	13 µL	50 µL
Rapamycin (Autophagy Inducer)	25 nmol	25 nmol
10X Assay Buffer	8 mL	30 mL
Chloroquine Control	7.5 µmol	7.5 µmol

ADDITIONAL MATERIALS REQUIRED

- Flow cytometer equipped with 488 nm laser source
- Standard fluorescence microscope
- Fluorescence microplate reader
- Tubes appropriate for holding cells for the flow cytometer
- Calibrated, adjustable precision pipettors, preferably with disposable plastic tips
- Adjustable speed centrifuge with swinging buckets (for suspension cultures)
- Deionized water
- Anhydrous DMSO
- Total growth medium suitable for cell type
- Indicator-free cell growth medium (e.g. Sigma Prod. No. D5921)
- EBSS (Sigma Prod. No. E2888)
- FBS (Fetal Bovine Serum)
- Glass microscope slides
- Glass cover slips of appropriate size
- 96-well tissue culture microplate with black wall and clear bottom

METHODS AND PROCEDURES

The procedures described in this manual assume that the user is familiar with the basic principles and practices of flow cytometry and is able to run samples according to the operator's manual pertaining to the instrument being used



Gently
hand mix
or vortex

NOTE: Allow all reagents to thaw at room temperature before starting with the procedures. Once thawed, gently hand-mix or vortex the reagents prior to use, to ensure a homogenous solution. Briefly centrifuge the vials at the time of first use, as well as for all subsequent uses, to gather the contents at the bottom of the tube.

A. Reagent Preparation

1. Positive Controls

a) Rapamycin

Rapamycin, a lipophilic macrolide antibiotic, is a widely used inhibitor of the mammalian target of Rapamycin (mTOR) and a well-established inducer of autophagy in a diverse range of cell lines from yeast to mammalian cells including neuron-like cells. Rapamycin-induced autophagy is characterized by the accumulation of autophagic vacuoles and the stimulation of autophagic flux.

Rapamycin included in the kit is supplied lyophilized (25 nmol). Before using, resuspend lyophilized Rapamycin in 50 μ L of DMSO. Resulting 500 μ M stock solution can be further diluted in cell culture medium to a desired concentration (recommended starting concentration of Rapamycin as a positive control is 500 nM). However, the optimal final concentration is cell-dependent and should be determined experimentally for each cell line being tested. Unused stock Rapamycin may be stored in small aliquots at -20°C for one month.

b) Chloroquine

Chloroquine is an anti-inflammatory drug used in the treatment or prevention of malaria. It suppresses inflammation by increasing lysosomal pH, thereby inhibiting lysosomal activity.

Chloroquine provided in the kit may be used in combination with Rapamycin or starvation in monitoring autophagic flux. It is supplied lyophilized (7.5 μ mol) and should be centrifuged briefly to gather the material at the bottom of the tube. Reconstitute the lyophilized material in 125 μ L deionized water for a 60 mM stock solution. Depending on the applications and specific cell lines, it is recommended that treatment with the agent will be performed using 10-120 μ M final concentration in order to observe changes in autophagic flux. Unused stock Chloroquine may be stored in small aliquots at -20°C for one month.

c) Rapamycin and Chloroquine as positive control

Some positive control cells (e.g., Jurkat cells) should be treated with Rapamycin and Chloroquine for 16-18 hours.

Recommended starting concentration of Rapamycin is 500 nM, Chloroquine is 10 μ M. Response to Rapamycin and Chloroquine is time and concentration dependent and may also vary significantly depending upon cell type and cell line. Negative control cells should be treated with a vehicle (DMSO, deionized water, media or other solvent used to reconstitute or dilute an inducer or inhibitor) for an equal length of time under similar conditions.

d) Starvation and Chloroquine as positive control

Some positive control cells (e.g., HeLa cells) should be starved in EBSS media in the presence of Chloroquine for 3-5 hours. Recommended starting concentration of Chloroquine is 10 μ M. Response to starvation and Chloroquine is time and concentration dependent and may also vary significantly depending upon cell type and cell line. Negative control cells should be incubated in the complete culture media with a vehicle (DMSO, deionized water, media or other solvent used to reconstitute or dilute Chloroquine) for an equal length of time under similar conditions.

2. 1X Assay Buffer

Allow the 10X Assay Buffer to warm to room temperature. Make sure that the reagent is free of any crystallization before dilution. Prepare enough 1X Assay Buffer for the number of samples to be assayed by diluting each milliliter (mL) of the 10X Assay Buffer with 9 mL of deionized water.

3. CYTO-ID® Green Detection Reagent 2

For optimal staining, the concentration of the CYTO-ID® Green dye 2 will vary depending upon the application. Suggestions are provided to use as guidelines, though some modifications may be required depending upon the particular cell type employed in the application.

a) **Fluorescence microscopy application:** Prepare a sufficient amount of **Microscopy Dual Detection Reagent** for the number of samples to be assayed as follows: For every 1 mL of 1X Assay Buffer or complete cell growth medium, add 2 μ L of CYTO-ID® Green Detection Reagent 2 and 1 μ L of Hoechst 33342 Nuclear Stain. If 1X Assay Buffer is used, supplement it with 5% FBS.

Note: *The dyes may be combined into one staining solution or each may be used separately, if desired.*



Important
Information

- *The Hoechst 33342 Nuclear Stain can be diluted further if its staining intensity is much stronger than that of the CYTO-ID® Green Detection Reagent 2. When staining BFP- or CFP-expressing cells, the Hoechst 33342 Nuclear Stain should be omitted due to its spectral overlap with these fluorescent proteins.*
 - *CYTO-ID® dye emits in the green region of the visible light spectrum and is thus not compatible with GFP.*
- b) **Flow cytometry application: Make a dye stain solution** by diluting 1 µL CYTO-ID® Green Detection Reagent 2 to 1 mL 1X Assay Buffer or cell culture medium without Phenol Red Indicator, supplemented with 5% FBS. For each sample to be stained, 250 µL of diluted CYTO-ID® Green Detection Reagent 2 staining solution will be used.
- c) **Fluorescence microplate application: Prepare a sufficient amount of Microplate Dual Detection Reagent** for the number of samples to be analyzed as follows: Add 1 µL of CYTO-ID® Green Detection Reagent 2 and 1 µL Hoechst 33342 Nuclear Stain into 1 mL 1X Assay Buffer or cell culture medium without Phenol Red Indicator, supplemented with 5% FBS. For each sample to be stained, 100 µL of Microplate Dual Detection Reagent will be used.

B. Live Cell Analysis by Fluorescence/Confocal Microscopy (Adherent Cells)

1. Grow cells on coverslips or tissue culture treated slides. When the cells have reached 50-70% level of confluence, carefully remove the medium.



Important
Information

NOTE: *Cells should be healthy and not overcrowded as results of the experiments will depend significantly on the cells' condition.*

2. Treat the cells with the testing reagent according to your experimental procedure. It is highly recommended to set up positive and negative controls within the same experiment (see *Section A, Step 1*).
3. Post-treatment, remove the medium with the testing reagents and positive control and wash the cells twice with 1X Assay buffer.

NOTE: *Be careful during washing procedure since autophagic cells can be easily dislodged from the slides. To preserve the cells, 2% - 5% FBS also may be added to the assay buffer at this point.*



Protect
from light



Microscopy

4. Dispense 100 μ L of Microscopy Dual Detection Reagent (see *Section A, Step 3*) to cover each sample of monolayer cells.
5. Protect samples from light and incubate for 30 minutes at 37°C.
6. Carefully wash the cells with 100 μ L of 1X Assay Buffer (see *Section A, Step 2*). Remove excess buffer and place coverslip on microscope slide.
7. An optional fixation step may be included at this step. Incubate for 20 minutes with 4% formaldehyde. Wash 3 times with 1X Assay Buffer.
8. Analyze the stained cells by wide-field fluorescence or confocal microscopy (60X magnification is recommended). Use a standard FITC filter set for imaging the autophagic signal. Optionally, image the nuclear signal using a DAPI filter set.

C. Live Cell Analysis by Fluorescence/Confocal Microscopy (Suspension Cells)

1. Cells should be cultured to a density not to exceed 1×10^6 cells/mL. Ensure that cells are in the log phase of growth before starting an experiment.



Important
Information

NOTE: *Cells should be healthy and not overcrowded since results of the experiments will depend significantly on the cells' overall condition.*

2. Collect the cells by centrifugation (5 minutes, 1000 rpm at room temperature). Resuspend the cells to a density of 1×10^6 cells/mL.
3. Treat the cells with the testing reagent according to your experimental procedure. It is highly recommended to set up positive and negative controls within the same experiment (see *Section A, Step 1*).
4. Post-treatment, remove the medium with the testing reagents and positive control and wash the cells twice with 1X Assay Buffer.
5. Carefully remove the supernatant and dispense 100 μ L of Microscopy Dual Detection Reagent solution (see *Section A, Step 3*) to cover the cell pellet. Resuspend the pellet by gently pipetting up and down.
6. Protect samples from light and incubate for 30 minutes at 37°C.
7. Wash the cells with 1X Assay Buffer. Remove excess buffer and re-suspend cells in 100 μ L 1X Assay Buffer.



Centrifuge

- An optional fixation step may be included at this step. Incubate for 20 minutes with 4% formaldehyde. Wash 3 times with 1X Assay Buffer.
- Apply a drop of the cell suspension onto a glass microscope slide and overlay with a coverslip.
- Analyze the stained cells by wide-field fluorescence or confocal microscopy (60X magnification recommended). Use a standard FITC filter set for imaging the autophagic signal. Image the nucleus using a DAPI filter set (optional).



Microscopy

D. Live Cell Analysis by Flow Cytometry

- Cells should be maintained via standard tissue culture practice. Grow cells overnight to log phase in a humidified incubator at 37°C, 5% CO₂.

NOTE: Cells should be healthy and not overcrowded since results of the experiments will depend significantly on the cells' overall condition. Cell density should not exceed 1×10^6 cells/mL.

- Treat cells with compound of interest according to experimental protocol. Prepare negative control cells using vehicle treatment. It is highly recommended to set up positive and negative controls within the same experiment (see Section A, Step 1)
- At the end of the treatment, trypsinize (adherent cells), or collect cells by centrifugation (suspension cells). Samples should contain 1×10^5 to 1×10^6 cells/mL.
- Centrifuge at 1000 rpm for 5 minutes to pellet the cells. Wash the cells by re-suspending the cell pellet in cell culture medium, 1X Assay Buffer, or other buffer of choice and collect the cells by centrifugation.
- Resuspend each live cell sample in 250 μ L of 1X Assay Buffer or indicator free cell culture medium containing 5% FBS.
- Add 250 μ L of the diluted CYTO-ID[®] Green stain solution (see Section A, Step 3) to each sample and mix well. Incubate for 30 minutes at room temperature or 37°C in the dark. It is important to achieve a mono-disperse cell suspension at this step by gently pipetting up and down repeatedly.
- After treatment, collect the cells by centrifugation and wash with 1X Assay Buffer. Resuspend the cell pellets in 500 μ L of fresh 1X Assay Buffer.
- An optional fixation step may be included at this step. Incubate for 20 minutes with 4% formaldehyde (or 10% formalin). Wash



Important Information

3 times with 1X Assay Buffer.

9. Analyze the samples in green (FL1) or orange (FL2) channel of a flow cytometer.

E. Live Cell Analysis by Fluorescence Microplate Reader

The CYTO-ID® Autophagy Detection Kit 2.0 has been shown to work for microplate readers. However, the conditions used for microscopy and flow cytometry may require additional optimization depending on cell line and end user applications.

For adherent cells

The procedure described below was developed using HepG2 and HeLa cells for which it is recommended that cells be seeded on plates at a density of 2.5×10^5 to 3.0×10^5 cells/mL, using 100 μ L cells/well. Any cell number and plate coating requirements should be optimized for the chosen cell model.

1. Seed cells in 96-well microplates, using 100 μ L cells/well, the day before the experiment, and allow cells to attach overnight under standard tissue culture practices. Cells should reach about 90% confluency to form a uniform monolayer in the well at the end of the experiment.

NOTE: *Cells should be healthy and not overcrowded as results of the experiments will depend significantly on the cells' overall condition.*

2. After overnight incubation, treat cells with compound of interest according to experimental protocol. Prepare negative control cells using vehicle treatment. It is highly recommended to set up positive and negative controls within the same experiment (see Section A, Step 1).
3. After the treatment, carefully remove the medium and dispense 100 μ L of 1X Assay Buffer to each well.

NOTE: *Be careful during washing procedure since autophagic cells can be easily dislodged from the plate. To preserve the cells, 2% - 5% FBS also may be added to the assay buffer at this point.*

4. Carefully remove all the buffer and dispense 100 μ L of Microplate Dual Detection Reagent (see Section A, Step 3) to each well.
5. Protect the sample from light and incubate at 37°C for 30 minutes.

6. Wash cells twice with 200 μ L of 1X Assay Buffer (see Note above) to remove excess dye and then add 100 μ L of 1X Assay Buffer to each well.
7. Analyze the plate with a fluorescence microplate reader. It is recommended to acquire data as soon after completing the assay as possible. The CYTO-ID[®] Green detection reagent 2 can be read with a FITC filter (Excitation \sim 480 nm, Emission \sim 530), and the Hoechst 33342 Nuclear Stain can be read with a DAPI filter set (Excitation \sim 340, Emission \sim 480). If the blue nuclear counterstain signal decreases by more than 30%, the compound is considered generally cytotoxic. Increases in the green autophagy signal after normalization with blue signal indicate the accumulation of the probe within the cells arising from an increase in autophagic vesicles.

For suspension cells

1. Culture the cells via standard tissue culture practice. Grow cells overnight to log phase in a humidified incubator at 37°C, 5% CO₂.
NOTE: *Cells should be healthy and not overcrowded since results of the experiments will depend significantly on the cells' overall condition. Cell density should not exceed 1x10⁶ cells/mL.*
2. Collect the cells by centrifugation (5 minutes, 1000 rpm at room temperature). Resuspend the cells to a density of 1x10⁶ cells/mL.
3. Treat cells with compound of interest according to experimental protocol. Prepare negative control cells using vehicle treatment. It is highly recommended to set up positive and negative controls within the same experiment (see Section A, Step 1).
4. At the end of the treatment, collect cells by centrifugation. Samples should contain 1x10⁵ to 1x10⁶ cells/mL. Wash the cells by re-suspending the cell pellet in cell culture medium, 1X Assay Buffer, or other buffer of choice and collect the cells by centrifugation. It is recommended that each suspension cell sample is tested in triplets. Following procedures are described accordingly.
5. Resuspend each cell sample in 400 μ L of CYTO-ID[®] Green Detection Reagent 2 solution (see Section A, Step 3). Incubate the cells for 30 minutes at 37°C in the dark. It is important to

achieve a mono-disperse cell suspension at this step by gently pipetting up and down repeatedly.

6. Wash the cells with 1X Assay Buffer. Remove excess buffer and re-suspend cells in 1X Assay Buffer. Count the cells and adjust the cells to a density of 5×10^5 cells/mL. If the number of the cells with testing reagent decreases by more than 30% compared to control, the compound is considered generally cytotoxic.
7. Add 100 μ L/well of the above cell suspension (e.g., 5×10^4 cells/well) to a 96-well microplate in triplicate, and analyze the cells with a fluorescence microplate reader. It is recommended to acquire data as soon after completing the assay as possible. The stain can be read with a FITC filter (Excitation ~480 nm, Emission ~530nm). Nuclear counterstain with Hoechst 33342 is optional for suspension cells, as the cell number has been normalized before adding to each well. Increases in the green autophagy signal indicate the accumulation of the probe within the cells arising from an increase in autophagic vesicles.

APPENDICES

Fluorescence Channel Selection

The selection of optimal filter sets for a fluorescence microscopy application requires matching the optical filter specifications to the spectral characteristics of the dyes employed in the analysis (see **Figure 2**). Consult the microscope or filter set manufacturer for assistance in selecting optimal filter sets for your microscope.

For flow cytometry, fluorescence channel FL1 (green) or FL2 (orange) is recommended for analysis of the CYTO-ID® Green Detection Reagent 2 staining using a 488nm laser source.




Figure 2. Absorbance and fluorescence emission spectra (499/548 nm) for CYTO-ID® Green Detection Reagent 2 (panel A) were determined in PBS. Absorbance and fluorescence emission spectra (350/461 nm) for Hoechst 33342 (panel B) were determined in 1X Assay Buffer.

Typical Results

A number of methods have been devised to investigate the autophagy pathway and the steps involved in the maturation of autophagosomes to autolysosomes, acid hydrolase-rich organelles in which the sequestered cytoplasmic material is ultimately degraded.⁽⁴⁻⁶⁾ For example, monodansylcadaverine (MDC) has been determined to be a useful probe for the analysis of the autophagic process by fluorescence microscopy.⁽⁵⁾ However, this probe requires 365nm UV illumination and thus, is not compatible with 488nm excitation sources commonly implemented in flow cytometry. The CYTO-ID® Autophagy Detection Kit 2.0 employs a 488nm-excitable green-emitting fluorescent probe to highlight the various vacuolar components of the autophagy pathway.

It should be noted that unlike the lysomotrophic dyes, MDC, LysoTracker® Red and Acridine Orange, which primarily detect lysosomes, the CYTO-ID® Green Detection Reagent 2 only weakly stains lysosomes, while serving both as a selective marker of autolysosomes and earlier autophagic compartments. This staining pattern differs markedly from that achieved with Lyso-ID® Red dye as well, which detects autophagolysosomes generated by Chloroquine and bafilomycin A treatment, but not vacuoles associated with other stimuli, such as serum starvation.⁽⁷⁾

Under physiological conditions, autophagy is a constitutive self-degradative process involved both in basal turnover of cellular components and as an induced response to nutrient starvation in eukaryotes. During autophagy, portions of the cytoplasm are sequestered by elongation of double-membrane structures called phagophores, which form vesicles called autophagosomes. These vesicles then fuse with lysosomes to form autolysosomes, where their contents are degraded by acidic lysosomal hydrolases for subsequent recycling (**Figure 1**). A prominent mammalian protein known to specifically associates with the autophagosome membrane is LC3-II.

When CYTO-ID® Green Detection Reagent 2 is incorporated into cells, the accumulation of this fluorescent probe is typically observed in spherical vacuoles in the perinuclear region of the cell, in foci distributed throughout the cytoplasm, or in both locations, depending upon the cell type under investigation.

In addition, an enhancement in the fluorescence emission intensity of the CYTO-ID® dye occurs upon compartmentalization with the lamellar membrane structures associated with autophagic vesicles. Oeste, et al., have reported the dye's selectivity for autophagic vesicles with BAECs transduced with RFP-LC3.

Besides Rapamycin treatment, there are several other approaches known to induce autophagy. One of the most potent known physiological inducers of autophagy is starvation. Autophagy induction can be observed with the CYTO-ID® Green Detection Reagent 2 within 1 hour of serum removal in both the HepG2 and HeLa cell lines. Another approach to activate autophagy is through the modulation of nutrient-sensing signal pathways. Several mTOR-independent autophagy activators have also been validated using the CYTO-ID® Autophagy Detection Kit 2.0 (**Table 1**). Lithium induces autophagy through inhibition of inositol mono-phosphatase (an mTOR-independent pathway). Trehalose and small-molecule enhancers of Rapamycin (SMERs) also induce autophagy by mechanisms that are not well

understood. Two FDA-approved compounds that induce autophagy in an mTOR-independent manner, Loperamide hydrochloride and Clonidine, also substantially increase green fluorescent signal in the assay.

Fluorescence Microscopy

Typical results of microscopy-based autophagy detection using the CYTO-ID[®] Green Detection Reagent 2 are presented in **Figure 3**.




Figure 3. HeLa cells were stained with CYTO-ID[®] Green Detection Reagent 2 after being cultured in (A) full media or (B) starvation media (EBSS) with 40 μ M Chloroquine for 4 hours. Cells starved in EBSS in the presence of Chloroquine showed very bright green fluorescent signals and punctate structures.

Typical results of flow cytometry-based analysis of cell populations using the CYTO-ID® Autophagy Detection Kit 2.0 are presented in **Figure 4**. In some cells lines (e.g., Jurkat cells), the formation and rapid degradation of autophagic vesicles by lysosomes may result in a low signal. In these cases, induction of autophagic flux is best visualized by accumulation of autophagic vesicles using rapamycin and lysosomal function is inhibited using chloroquine, preventing removal of these vesicles.




Figure 4. Flow cytometry-based profiling of autophagy in Jurkat cells. Jurkat cells were untreated or treated with 0.5 μ M Rapamycin (RAP), 10 μ M Chloroquine (CLQ) or both for 20 hours. After staining with CYTO-ID® Green Detection Reagent 2 for 30 minutes, cells were washed and analyzed by flow cytometry. Results are presented as histogram overlay. Cells treated with RAP + CLQ show an increase in fluorescence.

Fluorescence Microplate

Typical results of microplate-based autophagy detection using CYTO-ID® Green Detection Reagent 2 are presented in **Figure 5**. In some cells lines (e.g., Jurkat cells), the formation and rapid degradation of autophagic vesicles by lysosomes may result in a low signal. In these cases, induction of autophagic flux is best visualized by accumulation of autophagic vesicles using rapamycin and lysosomal function is inhibited using chloroquine, preventing removal of these vesicles.




Figure 5. Microplate-based profiling of autophagy in HepG2 cells. HepG2 cells were stained CYTO-ID® Green Detection Reagent 2 after being cultured for 20 hours in DMSO (control), 0.5 μ M Rapamycin (Rap), 10 μ M Chloroquine (CLQ), or both 0.5 μ M Rap and 10 μ M CLQ. Cells were also stained with Hoechst 33342 for cell number normalization. Cells treated with both Rapamycin and Chloroquine had an increase in autophagy, as measured by the CYTO-ID® Green Detection Reagent 2.

Autophagy Treatments

Treatment	Target	Effect	μM used	Induction Time (hr)	Cell Line
Starvation	Inhibits mammalian target of Rapamycin (mTOR)	Activates autophagy	N/A	1-4	HeLa, HepG2, Jurkat
Rapamycin	Inhibits mammalian target of Rapamycin (mTOR)	Activates autophagy	0.2	6-18	HeLa, Jurkat
PP242	ATP-competitive inhibitor of mTOR	Activates autophagy	1	18	HeLa
Lithium	Inhibits IMPase and reduce inositol and IP ₃ levels; mTOR-independent	Activates autophagy	10,000	18	HeLa, Jurkat
Trehalose	Unknown, mTOR-independent	Activates autophagy	50,000	6	HeLa, Jurkat
Bafilomycin A1	Inhibits Vacuolar-ATPase	Inhibits lysosome function	$6-9 \times 10^{-3}$	18	HeLa, Jurkat
Chloroquine	Alkalinizes Lysosomal pH	Inhibits lysosome function	10-50	18	HeLa, Jurkat
Tamoxifen	Increases the intracellular level of ceramide and abolishes the inhibitory effect of PI3K	Activates autophagy	4-10	6-18	HeLa, HepG2, Jurkat
Verapamil	Ca ²⁺ channel blocker; reduces intracytosolic Ca ²⁺ levels; mTOR-independent	Activates autophagy	40-100	18	HeLa, Jurkat
HydroxyChloroquine	Alkalinizes Lysosomal pH	Inhibits lysosome function	10	18	HeLa, Jurkat
Loperamide	Ca ²⁺ channel blocker; reduces intra-cytosolic Ca ²⁺ levels; mTOR-independent	Activates autophagy	5	18	HeLa
Clonidine	Imidazoline-1 receptor agonist; reduces cAMP levels; mTOR-independent	Activates autophagy	100	18	HeLa
MG-132	Selective proteasome inhibitor	Activates autophagy	2-5	18	HeLa, Jurkat
Norclomipramine	Alkalinizes Lysosomal pH	Inhibits lysosome function	5-20	18	HeLa
Epoxomicin	Selective proteasome inhibitor	Induce aggresome	0.5	18	HeLa
Velcade®	Selective proteasome inhibitor	Induce aggresome	0.5	18	HeLa
Amyloid beta peptide 1-42	Induce oxidative stress	Induce aggresome	25	18	SK-N-SH

Table 1. Treatments that influence autophagy, previously validated with CYTO-ID® Autophagy Detection Kit (ENZ-51031).

REFERENCES

1. Barth S, Glick D, Macleod KF. (2010) Autophagy: assays and artifacts. *J Pathol.* 221: 117-124.
2. Iwai-Kanai E, et al. (2008) A method to measure cardiac autophagic flux in vivo. *Autophagy.* 4: 322-329.
3. Chan L, et al. (2012) A novel image-based cytometry method for autophagy detection in living cells. *Autophagy.* 8: 1371-1382.
4. Mizushima, Yoshimori and Levine (2010) Methods in Mammalian Autophagy Research. *Cell.* 140: 313-326.
5. Vázquez and Colombo (2009) Assays to Assess Autophagy Induction and Fusion of Autophagic Vacuoles with a Degradative Compartment, Using Monodansylcadaverine (MDC) and DQ-BSA. *Methods in Enzymology, Volume 452:* 85-95.
6. Klionsky *et al* (2008) Guidelines for the use and interpretation of assays for monitoring autophagy in higher eukaryotes. *Autophagy* 4:2, 151-175.
7. Coleman *et al.* (2010) A live-cell fluorescence microplate assay suitable for monitoring vacuolation arising from drug or toxic agent treatment. *J Biomol Screen.* 15(4):398-405.
8. Oeste, *et al* (2013) Interactions between autophagic and endo-lysosomal markers in endothelial cells. *Histochem. Cell. Biol.* 139: 659-670.
9. Rossi *et al.* (2009) Desmethylclomipramine induces the accumulation of autophagy markers by blocking autophagic flux. *J Cell Sci* 122: 3330-3339.
10. Held *et. al.* (2011) A fluorescence microplate-based assay workflow for high-throughput screening of compounds modulating autophagy in living cells. BioTek Instruments, Inc. application note.
11. Kauntz H, Bousserouel S, Gossé F, Raul F. (2011) Silibinin triggers apoptotic signaling pathways and autophagic survival response in human colon adenocarcinoma cells and their derived metastatic cells. *Apoptosis.* 16:1042-53.
12. Klappan AK, Hones S, Mylonas I, Brüning A. (2012) Proteasome inhibition by quercetin triggers macroautophagy and blocks mTOR activity. *Histochem Cell Biol.* 137:25-36.
13. Lee JS, Lee GM. (2012) Monitoring of autophagy in Chinese hamster ovary cells using flow cytometry. *Methods.* 56: 375–382.



Product Manual

TROUBLESHOOTING GUIDE

Problem	Potential Cause	Suggestion
Low CYTO-ID® Green Detection Reagent 2 staining in all treatments, including positive control.	CYTO-ID® Green Detection Reagent 2 is photobleaching	Use mounting medium that prevents photobleaching. Optimize handling of the samples for fluorescence microscopy.
	A low concentration of the CYTO-ID® Green Detection Reagent 2 was used	Increase the reagent concentration (500X dilution of the dye is recommended)
	The incubation time with the dye reagent was insufficient	Increase the incubation time
	Concentration and/or time of treatment with autophagy inducer(s) is not optimized	The optimal final concentration of autophagy inducers (including positive control Rapamycin) is cell-dependent and should be determined experimentally for each cell line being tested.
High CYTO-ID® Green Detection Reagent 2 staining observed in negative control sample.	Cell cultures overgrown.	Suspension cells should not exceed a density 1×10^6 cells/mL and adherent cells should be approximately 50 - 70% confluent.
	Cell culture medium was depleted of nutrients.	Change media 4-8 hours before the experiment.
	Pathogen infection (Mycoplasma, etc.).	Obtain fresh cultures from reputable cell repository.

Problem	Potential Cause	Suggestion
The number of CYTO-ID® Green Detection Reagent 2 stained cells in the sample is too low to be readily quantified after assay.	Cell density/number was too low in the sample before the assay	Increase density/number of the cells in the sample
	Majority of the cells were lost during assay.	Autophagic cells may be loosely attached, so all staining and washing procedures should be performed gently. FBS or BSA (2-10%) may be added to the assay buffer.
CYTO-ID® Green Detection Reagent 2 staining fails to stain fixed and/or permeabilized cells.	The dye is only suitable for live-cell staining.	Use the dye for live-cell staining only. Cells can be fixed post staining as described.
Precipitate is observed in the 10X Assay Buffer	Precipitate forms at low temperatures.	Allow solution to warm to room temperature or 37°C, then vortex to dissolve all precipitate.
Cells do not appear healthy by microscopic examination post-assay.	Some cells require serum to remain healthy.	Add serum of BSA (2-10%) to the detection reagent and wash solutions.
Positive control (Rapamycin-treated) cells appear to be dead or are no longer attached to the plate surface.	The concentration and/or time of treatment are not optimized	The optimal final concentration of positive control (Rapamycin) is cell-dependent and should be determined experimentally for each cell line being tested.



NOTES

Product Manual



Product Manual

GLOBAL HEADQUARTERS

Enzo Life Sciences Inc.
10 Executive Boulevard
Farmingdale, NY 11735
Toll-Free:1.800.942.0430
Phone:631.694.7070
Fax: 631.694.7501
info-usa@enzolifesciences.com

EUROPE

Enzo Life Sciences (ELS) AG
Industriestrasse 17
CH-4415 Lausen
Switzerland
Phone:+41/0 61 926 89 89
Fax:+41/0 61 926 89 79
info-ch@enzolifesciences.com

For local distributors and detailed product information visit us online:

www.enzolifesciences.com