

TAISHAN泰山®


骨增量超级套装

坚如磐石 稳如泰山



骨钉 (金属锁定接骨螺钉)

用于固定骨块或骨片




编码	L
3301060215006	6
3301060215008	8
3301060215010	10
3301060215012	12
3301060215014	14
3301060215016	16

注: 8-14mm骨钉为套装内标配
另有6mm和16mm骨钉可单独采购

骨钻

制备帐篷钉和骨钉窝洞及皮质骨滋养孔



编码	Φd	L	L1
CCA.00017	Φ1.2	30	16
CCA.00018	Φ1.2	38	16
CCA.00019	Φ0.8	30	6.5
CCA.00020	Φ0.8	38	6.5

取骨钻

取环状骨块及受植区骨床预备



编码	L	Φd	ΦD	Φd1
CCA.00024	27.5	5.0	4.0	1.8
CCA.00025	27.5	6.0	5.0	1.8
CCA.00026	27.5	7.0	6.0	1.8
CCA.00027	27.5	8.0	7.0	1.8

膜钉 (平头固位钉)

用于固定屏障膜




L4

编码	L
CCA.00012	4

扁头帐篷钉 (伞状固位钉)

用于成骨空间维持



编码	L
CCA.00009	8
CCA.00010	10
CCA.00011	12

圆头帐篷钉 (伞状固位钉)

用于成骨空间维持



编码	L
CCA.00006	8
CCA.00007	10
CCA.00008	12

机用十字扳手

用于帐篷钉和骨钉的植入和取出



编码	L
CCA.00022	28
CCA.00023	35

手用十字扳手


用于帐篷钉和骨钉的植入和取出



编码	L
CCA.00021	80

快换手柄

用于连接手用十字扳手



编码	ΦD	L
CCA.0013	Φ3	103

起钉器


一侧用于取出膜钉，另一侧用于取出骨环



编码	L
CCA.00014	189

持钉器-弯


用于膜钉植入



编码	L
CCA.00015	181.5

持钉器-直

用于膜钉植入



编码	L
CCA.00016	150.5



微信公众号



微信订阅号



山东恒泰医疗器械有限公司
Shandong HENGTIME Medical Equipment Co., Ltd



www.hengtime.net



0538-8161828



山东省泰安市岱岳区南留街15号